

Project Name: Agobi Animation Page: 1 of 2



1-1 **REVISED (add front port and C02 caps to the product)** Product Spins, Text animated Transition-fade out


Text Box- The AGobiTM All-In-One Airway Stabilizing Device is a unique device featuring a revolutionary multifunctional adaptor providing ultimat stabilization for all endoscopic procedures, stabilization of airway patency in emergency situations, as well as other medical procedures requiring airway management.




2-3 Show Straps being put on

Text Box-Optional Adjustable Head Strap Included


Attach strap on one side, wrap around back of head to opposite side & secure



2-1 Transition-fade-in to head, zoom to mouth spin head slightlyback into zoom




2-3 show mouth opening and device being inserted into mouth & depressing tongue down stopping just short of vocal cords.



2-4 Show Connecting to High Flow Oxygen Source*
Text Box- Large Side port for Connect to High Flow
Oxygen Source*

*Multi-fuctional Adaptor can also be connected to Ambubags or Noninvasive Mechanical Ventilators Devices




2-5 **REVISED** (Show Co2 cap coming off) and tube being inserted for End Tidal CarbonDioxide Levels Monitoring

Text Box-Connect ETCO2 Monitoring Tube to small port




Project Name: Agobi Animation Page: 2 of 2



2-8 **REVISED** (Show front port cap opening) and endoscopy tube start to be inserted into front of device.


Text Box- Insert endoscope through front port

Front port can also be used for Suctioning Airway Secretions During or Post Sedation, or Anesthesia



2-9 show device extends to vocal cords- from inside the device

Text Box -End of Device extends to just short of vocal cords and esophagus for full tongue retraction and secure airway patency.



2-10 **REVISED** Switch views, show endoscope going down to the end of Agobi device, pause, (text transitions in/out, box out), continue animation of endoscope passing vocal cords into esophagus, fade-out. (with light on)

Text Box- The front port and tube's sleek design allows insertion of endoscope with ease.



1-1 Product Image Text animated in Transition-fade white



